



Simple Solution for Your Research

Hito Mini NisslStain™ Kit

[Catalog Number: HTKNS1020]

An easy to use Nissl staining system for the morphological characterization of the neurons and glial cells

User Manual And Material Safety Data Sheet

FOR IN VITRO RESEARCH USE ONLY

Hitobiotec Corp.

Simple solution for your research

Hito Mini NisslStain™ Kit

[Catalog Number: HTKNS1020]

An easy to use Nissl staining system for the
morphological characterization of the
neurons and glial cells

**User Manual
And
Material Safety Data Sheet**

FOR IN VITRO RESEARCH USE ONLY

Hitobiotec Corp.

© 2017 All Rights Reserved

Index

I.	Introduction	2
II.	Kit Contents	3
III.	Tissue Preparation	4
IV.	Staining Procedure	8
V.	Material Safety Data Sheet (MSDS)	9

I. Introduction

Hito mini NisslStain™ Kit is formulated for the staining of both neurons and glial cells. This solution can be used with frozen or paraffin-embedded tissue sections fixed with any fixative (formalin preferred).

Hito Mini NisslStain™ Kit is designed in a ready-to-use format and offers high quality, rapid and reproducible staining. Our kit is suitable for small scale application (for 40 or more slides).

Hito Mini NisslStain™ Kit has been tested extensively on the tissues from several species of animals and it is a simple solution for your research.

For photo samples, please visit our web site at
www.hitobiotec.com

II. Kit Contents

Store Hito Mini NisslStain™ Kit at room temperature

Stain Solution	50 ml
Differentiation Solution	50 ml
Staining Jar for Staining / differentiation(12 ml)	2
Dropping Bottle for Stain Solution	1
User Manual and MSDS	1



Note

Before using Hito Mini NisslStain™ Kit, please make sure you have the following **Required Equipment / Materials** in your lab (not included in the kit):

1. Cryostat or Microtome
2. Ethanol, xylene, double distilled or deionized water
3. Slide (recommend Hito Super-Safe Slide #HTHS0101)
4. Coverslips
5. Staining jars for slides wash
6. Resinous mounting medium
7. Light microscope

III. Tissue Preparation

For Frozen Section

1. Prepare perfusion system.
2. Prepare animal for infusion by administering a lethal dose of anesthesia. Monitor it until the point when the animal fails to respond to pinching of the foot.
3. Cut the skin of the mouse from the abdomen to the top of the thorax. Open the abdominal wall below the ribcage. Lift the sternum with tweezers and cut the diaphragm. Then cut away the lower part of the ribcage to partially expose the heart.
4. Quickly insert needle of infusion set into left ventricle. Clamp the needle in place.
5. Begin perfusion of PBS very slowly (i.e., 10 to 15 ml/min). After the perfusion system begins pumping the PBS, immediately cut the inferior vena cava to allow an escape route for the blood and perfusion fluid.
6. Perfuse PBS at a moderate to rapid rate (<20 ml/min) and continue until the effluent runs clear, which may require 25 to 50 ml of solution.
7. After the effluent runs clear, stop the pump and introduce 4% PFA into the infusion set line running into the animal. Perfuse 50-100 ml 4% PFA at a moderate to slow rate.
8. Remove and transfer tissue into 4% PFA, store at 4°C. Replace 4% PFA after 24 hours, and continue to store at 4°C for 24 hours.
9. Transfer the tissue into 20% sucrose solution, store at 4°C. Replace 20% sucrose solution after 24 hours, and continue to store at 4°C for 24 hours, until the tissue sinks into the sucrose solution.
10. Place ~300 to 500 ml isopentane in a metal container large enough to hold a corresponding sieve-like basket. Place the metal container with the isopentane in dry ice for 15 to 30 min, until the temperature of the isopentane reaches -70°C.

11. Place the tissue briefly on absorbent paper to remove excess solution.
12. Place the tissue in OCT compound in a base mold on the mesh bottom of the sieve-like basket.
13. Slowly immerse the basket with the tissue in the cooled isopentane for 30 sec. to 1 min.



Note

The time of immersion is absolutely critical; it must be long enough to result in complete freezing of the tissue, but not so long that the tissue cracks. It may be necessary to test various times to determine the one that is optimal to meet these criteria.

14. Wrap the dried, frozen tissue block in aluminum foil and store at -70°C until sectioning is performed.
15. Set the cryostat chamber temperature at -17°C .



Note

The -17°C setting is satisfactory in most cases, but may need optimization for different cryostat and tissue types in order to cut sections smoothly and keep integrity.

16. Place specimen holder / cryostat chuck on dry ice and add embedding matrix or water on the surface of the specimen holder / chuck. As the embedding matrix or water begins to freeze, place the frozen tissue block into it so that the tissue block adheres to the specimen holder / chuck.
17. Slowly cut the tissue into sections (11-25 μm thickness) on a cryostat with the chamber temperature set at -17°C .
18. Mount the sections on the Hito Super-Safe Slide
19. Air dry slides (30 minute) at room temperature. Dried sections should be processed as soon as possible but may be stored in a slide box at -20°C for one year.

For Paraffin Tissue Section

1. Prepare perfusion system.
2. Prepare animal for infusion by administering a lethal dose of anesthesia. Monitor it until the point when the animal fails to respond to pinching of the foot.
3. Cut the skin of the mouse from the abdomen to the top of the thorax. Open the abdominal wall below the ribcage. Lift the sternum with tweezers and cut the diaphragm. Then cut away the lower part of the ribcage to partially expose the heart.
4. Quickly insert needle of infusion set into left ventricle. Clamp the needle in place.
5. Begin perfusion of PBS very slowly (i.e., 10 to 15 ml/min). After the perfusion system begins pumping the PBS, immediately cut the inferior vena cava to allow an escape route for the blood and perfusion fluid.
6. Perfuse PBS at a moderate to rapid rate (<20 ml/min) and continue until the effluent runs clear, which may require 25 to 50 ml of solution.
7. After the effluent runs clear, stop the pump and introduce 4% PFA into the infusion set line running into the animal. Perfuse 50-100 ml 4% PFA at a moderate to slow rate.
8. Remove and transfer tissue into Hito Bouin's Plus solution, store at 4°C. Replace Hito Bouin's Plus solution after 24 hours, and continue to store at 4°C for 24-48 hours.
9. After fixation, dehydrate the tissue in a graded ethanol/water series at room temperature: (for 0.5 cm³ tissue)
 - 70% ethanol for 2 changes of 2 hours each
 - 95% ethanol for 2 changes of 1.5 hour each
 - 100% ethanol for 2 changes of 1 hour each
10. Replace ethanol with xylene for 2 changes, each 45 minutes at room temperature.

11. Immerse the tissue in the paraffin wax (56-58°C), 2 changes, 1.5 hour each.



Note

The time of immersion is absolutely critical; it must be long enough to remove xylene in the tissue, but not so long that the tissue stiffens. It may be necessary to test various times to determine the one that is optimal to meet these criteria.

12. Embed tissues in paraffin blocks
13. Turn on the water bath and check that the temperature is 45°C. Use fresh deionized water. Insert the block into the microtome chuck. Set the dial to cut 5-11 μm sections. Cut sections and pick them up with forceps or a fine paint brush and float them on the surface of the water bath. Float the sections onto the surface of Histo Super-Safe Slide.
14. Place the slides with paraffin sections in a 60°C oven for 2 hours (so the wax just starts to melt) to bond the tissue to the glass. Slides can be stored in slide box at room temperature.

IV. Staining Procedure

(For Frozen Section and Paraffin Tissue Section)

1. Place slides in xylene 2 times, 3 minutes each.
2. Place slides in 100% ethanol 2 times, 3 minutes each.
3. Place slides in 95% ethanol 2 times, 3 minutes each.
4. Place slides in 75% ethanol 2 times, 3 minutes each.
5. Place slides in 50% ethanol 2 times, 3 minutes each.
6. Rinse slides in distilled water for 1 minute.
7. Using a 12 ml staining jar, Place slides in Stain Solution, (or from dropping bottle place a few drops of Stain Solution on the section) and allow to stand for 0.5 - 5 minutes depending on the desired intensity.
8. Rinse slides in distilled water for 1 minute.



Note

if the tissue is overstained, the tissue needs to be differentiated with Differentiation Solution, 1-5 dips. Differentiation requires some practical experience to ascertain the correct end-point. The staining intensity of both cellular elements and background will fade rapidly in this solution. Check under the microscope and if over differentiated, return to step 6.

9. Rinse slides in distilled water for 1 minute.
10. Air dry tissue slides (1-3 hours or overnight) at room temperature
11. Clear in xylene, 2 times, 3 minutes each, and apply coverslip over sections using xylene based resinous mounting medium.
12. Allow to dry. The slide can be viewed after drying by bright field microscopy.

V. Material safety data sheet (MSDS)

Date Updated: 11/01/2016
Version 1.5

1. Product and Company Information

Product Name	Hito Mini NisslStain™ Kit
Product Number	HTKNS1020
Brand	Hitobiotec
Company Address	Hitobiotec Corp. P.O.Box 7528 Kingsport, TN 37664 USA
Technical Phone:	423-520-6880
Emergency Phone:	423-520-6880

2. Composition and Information on Ingredient

Substance Name	CAS #	SARA 313
Hito Mini NisslStain™ Kit	None	No

Ingredient Name	CAS #	SARA 313
WATER	7732-18-5	No
Sodium Acetate	6131-90-4	No
Cresyl Violet	10510-54-0	No
Acetic acid	64-19-7	Yes

3. Hazards Identification

EMERGENCY OVERVIEW

Harmful if swallowed. Causes skin irritation. Causes serious eye irritation. May cause respiratory irritation

HMIS RATING

HEALTH: 1 FLAMMABILITY: 0 REACTIVITY: 0

NFPA RATING

HEALTH: 1 FLAMMABILITY: 0 REACTIVITY: 0

Potential Health Effects

Inhalation May be harmful if inhaled. Causes respiratory tract irritation.

Skin May be harmful if absorbed through skin. Causes skin irritation.

Eyes Causes eye burns, eye irritation.

Ingestion Toxic if swallowed. May be fatal if swallowed. Causes burns.

4. FIRST AID MEASURES

General advice

Consult a physician. Show this safety data sheet to the doctor in attendance.
Move out of dangerous area.

If inhaled

If breathed in, move person into fresh air. If not breathing give artificial respiration

In case of skin contact

Take off contaminated clothing and shoes immediately. Wash off with soap and plenty of water. Consult a physician.

In case of eye contact

Continue rinsing eyes during transport to hospital. Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.

If swallowed

Do NOT induce vomiting. Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.

5. FIRE-FIGHTING MEASURES

Flammable properties

Flash point no data available

Ignition temperature no data available

Suitable extinguishing media

Use water spray, alcohol-resistant form, dry chemical or carbon dioxide.

Special protective equipment for fire-fighters

Wear self contained breathing apparatus for fire fighting if necessary.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions

Use personal protective equipment. Avoid dust formation. Avoid breathing dust. Ensure adequate ventilation. Evacuate personnel to safe areas.

Environmental precautions

Prevent further leakage or spillage if safe to do so. Do not let product enter drains. Discharge into the environment must be avoided.

Methods for cleaning up

Pick up and arrange disposal without creating dust. Keep in suitable, closed containers for disposal.

7. HANDLING AND STORAGE

Handling

Perform experiment in a properly functioning chemical hood, which is vented to the outside. Wear glasses and disposable gloves while handling kit reagents. Wash hands thoroughly after performing the test.

Storage

Keep container tightly closed in a dry and well-ventilated place. Store at room temperature, preferably in a cool place.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Contains no substances with occupational exposure limit values.

Personal protective equipment

Respiratory protection

Where risk assessment shows air-purifying respirators are appropriate use a full-face respirator with multipurpose combination (US) or type ABEK (EN 14387) respirator cartridges as a backup to engineering controls. If the respirator is the sole means of protection, use a full-face supplied air respirator. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Hand protection

Handle with gloves.

Eye protection

Safety glasses with side-shields conforming to EN166

Hygiene measures

Avoid contact with skin, eyes and clothing. Wash hands before breaks and immediately after handling the product.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance

Form liquid

Safety data

pH no data available
Melting point no data available
Boiling point no data available
Flash point no data available
Ignition temperature no data available
Lower explosion limit no data available
Upper explosion limit no data available
Water solubility no data available

10. STABILITY AND REACTIVITY

Storage stability

Stable under recommended storage conditions. Avoid Light.

Materials to avoid

Strong oxidizing agents, metals

Hazardous decomposition products

Hazardous decomposition products formed under fire conditions. - Hydrogen chloride gas, Mercury/mercury oxides, Potassium oxides, Chromium oxides.

11. TOXICOLOGICAL INFORMATION

Acute toxicity Refer to component MSDS

Irritation and corrosion Refer to component MSDS

Sensitisation Refer to component MSDS

Signs and Symptoms of Exposure

no data available

Potential Health Effects

Inhalation May be fatal if inhaled. Material is extremely destructive to the tissue of the mucous membranes and upper respiratory tract.

Skin Causes skin burns, skin irritation. May be fatal if absorbed through skin.

Eyes Causes eye burns, eye irritation.

Ingestion Toxic if swallowed. May be fatal if swallowed. Causes burns.

12. ECOLOGICAL INFORMATION

Elimination information (persistence and degradability)

Refer to component MSDS

Ecotoxicity effects

Refer to component MSDS

Further information on ecology

Refer to component MSDS

13. DISPOSAL CONSIDERATIONS

Product

Observe all federal, state, and local environmental regulations. Contact a licensed professional waste disposal service to dispose of this material.

Contaminated packaging

Dispose of as unused product.

14. TRANSPORT INFORMATION

DOT (US)

Not dangerous goods

IMDG

Not dangerous goods

IATA

Not dangerous goods

15. OTHER INFORMATION

Further information

The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. Hitobiotech, Inc., shall not be held liable for any damage resulting from handling or from contact with the above product. See Terms & Conditions page on our website for additional terms and conditions of sale.

Notes

Hitobiotec Corp.

P.O.Box 7528
Kingsport, Tennessee
U.S.A.

Phone: 423-520-6880

Email: info@hitobiotec.com

www.hitobiotec.com

